
Osteograft[®]

allogene transplantate

HERSTELLUNG & KOOPERATION:
PRODUCTION & COOPERATION:

PRODUKTKATALOG
2026/2027

INFORMATION & DISTRIBUTION:

■ **ARGON** Dental Vertriebs GmbH & Co. KG

Franz-Kirsten-Straße 1
D-55411 Bingen am Rhein

Fon +49 (0) 6721 / 3096-0
Fax +49 (0) 6721 / 3096-29

www.argon-dental.de
info@argon-dental.de

■ **DENTALES PROGRAMM**

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DENTAL PROGRAM

Product list Soft tissue management – Updated on 11/07/2023

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DENTAL TRANSPLANTS / SOFT TISSUE MANAGEMENT

Jaw and Maxillofacial Surgery /Stomatology

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E

DIZG – THE GERMAN INSTITUTE FOR CELL- AND TISSUE REPLACEMENT

Was founded as a non-profit institution in 1993 by experienced physicians and scientists from the long established tissue banks of the University of Leipzig and the Charité University Hospital in Berlin.

The Institute is dedicated to the issuance of tissue donations to the medical community. For patients requiring transplants in form of high quality and safe allograft tissue transplants.

Approximately 500,000 human hard and soft tissue transplants have been provided to hospitals and surgeons worldwide since its founding.

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E

In Germany, human cell and tissue transplants are regulated as pharmaceutical products (drugs) and require a manufacturing license and pharmaceutical approval by the national authorities and a marketing authorization from the competent authority respectively. The biological safety and all aspects of the manufacturing process as well as clinical application of the transplants are continually monitored and controlled by the relevant governmental authorities. DIZG has 10 manufacturing licenses for human tissue transplants.

Our biological safety system includes:

- An extensive donor screening
- A comprehensive serological testing scheme that exceeds the standards set by the
- EU-Directive 23/2004 and includes 4 Virus PCR tests
- A validated process for the deactivation and removal of viruses, bacteria and fungi
- In-process controls and testing of finished transplants

The DIZG and ARGON are DIN EN ISO 13485 certified and adhere to the ethical guidelines and quality standards of the European Association of Tissue Banks (EATB).

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GENERAL INFORMATION

E

The application of the grafts of the DIZG is subject to documentation according to the Transplantation Act (TPG). For each package, labels are provided for attachment to the patient's file.

All dimensions are approximate or minimum values. This is due to both the nature of the manufacturing process and inter-donor variance. Advice about the clinical application of the transplants in detail is found in the „instructions for use“ information provided in every package.

PRESERVATION METHODS

The transplants designated with a GT catalog number are freeze-dried. (Exception: Amnion transplants are dried.) Freeze-dried transplants have an expiry date of 5 years from the date of manufacture when stored in unopened packaging at or below 25° C. DBM Pasty is ready for immediate use without rehydration.

All other freeze-dried transplants must be rehydrated in a suitable physiological medium (e.g., isotonic saline) for at least 30 minutes.

*Advancing
the
World of
Tissue
Transplantation*

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D

APPLICATIONS / INDICATIONS

Class I to III, Hemisection, Furcation, Root Tip Resection, Extraction, Wisdom Tooth Removal, Sinus Elevation, Implantation and augmentation. The main application for allogenic tissue transplants is the defect filling in oral and maxillofacial surgery, especially in particular indications for which this type of material, to an ever-increasing extent, represents the only and most successful choice. DBM Putty and DBM X-press Mix can be used immediately due to the sodium hyaluronate. All other materials are to be rehydrated before use with a suitable sterile physiological medium (e.g., isotonic infusion). An ideal substrate for transplants is OsteoGel. OsteoGel is a plasma membrane protein carrier sodium hyaluronate. Sodium hyaluronate plays an essential role in cell proliferation, migration, and adhesion and provides very good osteogenesis. (Delivered in syringe to offer easy handling).

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E

DEMINERALIZED BONE MATRIX (DBM) IS:

- Fully demineralized human bone matrix
- Almost pure collagen (95%, 5% BMP)
- Radiopaque
- Approved as a drug according to AMG
- Corticocancellous and freeze-dried
- Available in two grain sizes and easily processable

HARTGEWEBEMANAGEMENT /
REKONSTRUKTIVE CHIRURGIE

HARD TISSUE MANAGEMENT / RECONSTRUCTIVE SURGERY

DEMINERALISIERTE KNOCHENMATRIX (DBM) PASTÖS*
DEMINERALIZED BONE MATRIX (DBM) PASTY

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E

DEMINERALIZED BONE MATRIX (DBM) OFFERS:

- Very fast remodeling times (< 3 months)
- Reliably vascularized bone without residues of foreign material
- Ready for immediate use for bone substitute demands (DBM Pastös/Putty)
- Wide range of applications for all minor defect fillings in oral surgery
- Very good usability within the Periodontology
- For enclosed defects

HARTGEWEBEMANAGEMENT / REKONSTRUKTIVE CHIRURGIE

HARD TISSUE / RECONSTRUCTIVE SURGERY

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E

CORTICAL, CANCELLOUS AND CORTICOCANCELLOUS GRANULES ARE

- Mineralized human bone matrix (10% mineral content)
 - average remodeling time from 6 to 8 months (corticospongiuous granules)
 - long remodeling time of 12 months (corticalis granules)
- Non-radiopaque
- Freeze-dried and approved as a drug according to AMG
- Consist of corticospongiuous bone or pure tubular bone
- Are available in two or three granulate sizes and easy to process

HARTGEWEBEMANAGEMENT / REKONSTRUKTIVE CHIRURGIE

HARD TISSUE / RECONSTRUCTIVE SURGERY

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E

CORTICAL, CANCELLOUS AND CORTICOCANCELLOUS GRANULES OFFER

- High resorption protection as maximum assurance for volume stability
- High predictable results for large augmentations
- The basis for membrane technology (Tent-Pole or Umbrella)
- The ideal filling material for augmentation in shell technique with fixed or flexible bone shields
- The ideal material for contouring of Bone Block Onlays

HARTGEWEBEMANAGEMENT /
REKONSTRUKTIVE CHIRURGIE

HARD TISSUE / RECONSTRUCTIVE SURGERY

DEMINERALISIERTE KNOCHENMATRIX (DBM) PASTÖS*
DEMINERALIZED BONE MATRIX (DBM) PASTY / GRANULES

DEMINERALISIERTE KNOCHENMATRIX (DBM) PASTÖS* / PUTTY
DEMINERALIZED BONE MATRIX (DBM) PUTTY / PASTY

(ohne Mineralien – Matrix für sofortige Zellen-Infiltration)
(without mineral - matrix for immediate cell infiltration)

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** Putty* is mixed with the plasma membrane protein carrier sodium hyaluronate. Sodium hyaluronate plays an essential role in cell proliferation, migration and adhesion and offers excellent osteogenesis. (Delivered in injection form to offer excellent handling).

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OSTEOGEL BARRIER, THE SAFE, SYNTHETIC AND USER-FRIENDLY ALTERNATIVE TO ABSORBABLE MEMBRANES.

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
E

After even application of the barrier gel and due to the relatively short 6- to 12-hour resorption time of the uncross-linked hyaluronic acid, a homogeneous layer of crosslinked particles remains on the surface of the defect. As the particles vary in size, they condense and form a flexible membrane on the area to be shielded. The barrier function is guaranteed for 3 weeks.

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FOR IMPROVED OR CONTROLLED WOUND HEALING
AFTER SURGICAL PROCEDURES.

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Osteogel Hyaluronic Acid Liquid Gel is locally absorbed through the oral mucosa within 6 to 12 hours. To achieve optimal treatment success, it is important that the hyaluronic acid liquid gel is applied to the wound and that it is then closed with a suture.

The low viscosity of the hyaluronic acid liquid gel enables easy clinical application.

D

E

HARTGEWEBEMANAGEMENT /
REKONSTRUKTIVE CHIRURGIE

HARD TISSUE / RECONSTRUCTIVE SURGERY

KNOCHEN BLOCK – KNOCHEN SPAN – KNOCHEN KEIL
BONE BLOCK – BONE PLATE – BONE WEDGE

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
E

APPLICATIONS / INDICATION

Block Graft Onlay and Container Techniques
Lateral augmentation with or without simultaneous implantation
Shell or block techniques

The main application of these allogenic hard tissue grafts is the massive bone augmentation in the course of major implantological operations as a reliable and patient-friendly alternative to autologous bone. Volume-stable and suitable for long remodeling times. All block, chip and wedge-shaped grafts should be rehydrated with a suitable sterile physiological medium (e.g. isotonic infusion solution).

Preferred filling materials are the granules corticalis and corticocancellous. All granules except DBM Putty should be rehydrated with a suitable sterile physiological medium (e.g. isotonic infusion solution).

An ideal carrier material to supplement the filling material is OsteoGel. OsteoGel is a plasma-glucose protein carrier sodium hyaluronate solution. Sodium hyaluronate plays an essential role in cell proliferation, migration and adhesion and provides ideal osteogenesis. (Supplied in injection form and offers optimal handling).

HARTGEWEBEMANAGEMENT /
REKONSTRUKTIVE CHIRURGIE

HARD TISSUE / RECONSTRUCTIVE SURGERY

CORTI-FLEX® / DEMINERALISIERTE KNOCHENMATRIX
CORTIFLEX® / DEMINERALIZED BONE MATRIX

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E

CORTI-FLEX IS

- Demineralized cortical chip of human origin
- Flexible after rehydration
- Drug approved according to AMG
- Freeze-dried and freeze-preserved
- Stable contour casing (thick variant)
- Bony membrane (thin variant)
- Resorption time thin 6 to 12 months
Resorption time thick 18 to 24 months


CORTI-FLEX OFFERS

- Flexible adaptation for contour shuttering
- Ideal for single-stage implantation and augmentation
- Good compression of the granulated material
- Easy fixation with titanium pins (thin) or with lag screws and washers

HARTGEWEBEMANAGEMENT / REKONSTRUKTIVE CHIRURGIE

HARD TISSUE / RECONSTRUCTIVE SURGERY

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E

FEMUR BONE PLATE IS

- Compact bone from the femur (thigh bone)
- Human origin and drug approved according to AMG
- Freeze-dried and freeze-preserved
- Long stable for resorption protection
- The ideal alternative to autografts
- Basis for the container technique with allogenic materials

FEMUR BONE PLATE OFFERS

- Stable shield for container technique with granulated material
- Safety for very long remodeling times
- Ideal conditions for lateral and vertical augmentation
- Sufficient thickness for individual manual adaption
- Reliable fixation with set screws

HARTGEWEBEMANAGEMENT /
REKONSTRUKTIVE CHIRURGIE

HARD TISSUE / RECONSTRUCTIVE SURGERY

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E

CANCELLOUS BONE BLOCKS, - CUBES, -
CYLINDERS ARE

- Human cancellous bone without cortical plate
- Freeze-dried and according to AMG pharmaceutically approved
- Thanks to the DIZG's own sterilization process very quickly revascularizable
- Easy to adapt
- Ideal for captured, non-exposed defects

CANCELLOUS BONE BLOCKS, - CUBES, -
CYLINDERS OFFER

- The base for the bone ring technique in Implantology
- An alternative defect filling for large Extraction sockets
- Inlays for the sinus
- Fast conversion times thanks to quick Revascularization
- Ideal primary stability with single-stage implantation

HARTGEWEBEMANAGEMENT / REKONSTRUKTIVE CHIRURGIE

HARD TISSUE / RECONSTRUCTIVE SURGERY

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E

CORTICOCANCELLOUS BONE BLOCKS ARE

- Human spongy blocks with one, two or three-sided cortical plates
- Originating from the os-ilium or femoral head
- Thanks to cortical shield and slow conversion reliably protected against resorption
- Ideal as an alternative to lifting of retromolar bone blocks or a painful invasive procedure on the iliac crest
- Freeze-dried
- Approved as a drug according to AMG

CORTICOCANCELLOUS BONE BLOCKS OFFER

- Predictable results for lateral and vertical augmentation
- Easy intraoperative adaptability
- The base for the Osteograph-3D-Technology
- Long remodeling times for gentle conversion with good revascularization
- A protective shield when dehiscence appears
- Various sizes and shapes for indication-based block selection

HARTGEWEBEMANAGEMENT/
REKONSTRUKTIVE CHIRURGIE

HARD TISSUE / RECONSTRUCTIVE SURGERY

CORTI-FLEX® / DEMINERALISIERTE KNOCHENMATRIX (DBM)
CORTI-FLEX® DEMINERALIZED FLEXIBLE CORTICAL BONE PLATE

FEMUR / KNOCHEN SPAN und KEIL
FEMUR / NON-FLEXIBLE CORTICAL BONE PLATE and CHUCK

CORTI-FLEX® / DEMINERALISIERTE KNOCHENMATRIX (DBM) SPAN
CORTI-FLEX® / DEMINERALIZED BONE MATRIX Flexible Bone Plate

(Konturverschalung ± 180 Tage)
(contour casing ± 180 days)

Artikelnummer Stock Number	Produktbeschreibung	Product Description	Größe Size
GT 3540	CortiFlex DBM Span	CortiFlex DBM flexible plate	15 x 50 x 2 mm
GT 3541	CortiFlex DBM Span	CortiFlex DBM flexible plate	15 x 100 x 2 mm
GT 3545	CortiFlex DBM Span	CortiFlex DBM flexible plate	15 x 30 x 2 mm
GT 3546	CortiFlex DBM Span	CortiFlex DBM flexible plate	15 x 60 x 2 mm
GT 3547	CortiFlex DBM Span	CortiFlex DBM flexible plate	15 x 120 x 2 mm
GT 3548	CortiFlex DBM Span	CortiFlex DBM flexible plate	20 x 25 x 2 mm
CORTI-FLEX® Dünn			
GT 3545T	CortiFlex dünn (ca. 0,5 mm)	CortiFlex thin (ca. 0.5 mm)	15 x 30 x 0,5 mm
GT 3548T	CortiFlex dünn (ca. 0,5 mm)	CortiFlex thin (ca. 0.5 mm)	20 x 25 x 0,5 mm

FEMUR / KNOCHEN SPAN und KEIL
FEMUR / BONE PLATE and CHUCK

(Schalentechnik ± 180 Tage)
(Shell technique ± 180 days)

GT 1601	Femur Span	Femur Plate	15 x 20 x 2 mm
GT 1603	Femur Span	Femur Plate	10-15 x 50 x 2 mm
GT 1606	Femur Span	Femur Plate	10 x 20 x 1 mm
GT 1607	Femur Span	Femur Plate	15 x 20 x 1 mm

HARTGEWEBEMANAGEMENT/
REKONSTRUKTIVE CHIRURGIE

HARD TISSUE / RECONSTRUCTIVE SURGERY

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E

APPLICATIONS / INDICATIONS

Barrier Membrane and soft tissue regeneration (GTR)

The main application of these allogenic tissue transplants is the defect filling in oral and maxillofacial surgery in special indications for which this type of material, to an ever increasing extent, represents the only and most successful material of choice.

Epiflex as well as fascia lata are to be prepared before use with a suitable sterile physiological medium (eg isotonic infusion solution) for rehydration.

An ideal carrier for regeneration is OsteoGel. OsteoGel is a plasma membrane protein carrier sodium hyaluronate solution. Sodium hyaluronate plays an essential role in cell proliferation, migration and adhesion and provides excellent osteogenesis. (Will be delivered in injection mold and offers excellent handling). The Osteogel Barrier is also used as a „liquid membrane“ and provides a long-term effect for soft tissue regeneration.

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E

FASCIA LATA IS

- Thick connective tissue cover (fascia)
- Human origin
- Comes from the thigh
- Freeze-dried
- Approved as a drug according to AMG
- Highly tear-resistant, stable and easy to adapt

FASCIA LATA OFFERS

- Exceptionally long remodeling time (> 6 months)
- Reliable protection of the graft
- Insensitive to suture dehiscence
- The base for Tent-Pole and Umbrella techniques
- Good and very stable intraoperative handling

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E

APPLICATIONS / INDICATIONS

Membrane and soft tissue regeneration (GTR) and also as a connective tissue transplant. The main application of these materials is as a membrane in oral and maxillofacial surgery, in particular in special indications for which this type of material, to an ever-increasing extent, represents a successful choice. The membrane is to be rehydrated before use with a suitable sterile physiological medium (eg isotonic infusion solution). An ideal carrier material for wound healing is OsteoGel. The OsteoGel is a plasma membrane protein carrier sodium hyaluronate solution. Sodium hyaluronate plays an essential role in cell proliferation, migration and adhesion and provides excellent osteogenesis. (Will be delivered in injection mold and offers excellent handling).

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EPIFLEX IS:

- Dermis of human origin
- Completely free of cells
- Very flexible and very tear-resistant
- Permeable to water vapor and oxygen
- Pharmaceutical approval
- Under the german medical legislation

EPIFLEX OFFERS:

- An ideal structure for resettlement by cells and vessels
- Sterility and excellent intraoperative handling
- Long shelf life at room temperature
- Excellent cosmetic results
- Very high compatibility

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E

Epiflex is the only acellular dermis of human origin that is pharmaceutically approved*

- Implant surgery
- Additional tissue strengthening
- For thin skin layers
- Suitable for one-time reconstruction
- Complication surgery
- Breast reconstruction after mastectomy

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E

Tensile strength (Tension) 70N
Strength at rupture 53N
Modulus of elasticity 15 MPa

BIOMETRICS OF EPIFLEX IS COMPARABLE TO NATIVE DERMIS

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Wasserdampf 1,7 mol m^{-2s-1}
Sauerstoff 2,7 mol m^{-2s-1}

Water vapor 1.7 mol m^{-2s-1}
Oxygen 2.7 mol m^{-2s-1}

THE PERMEABILITY PROPERTIES ARE COMPARABLE TO THOSE OF NATIVE DERMIS:

* Epiflex A new decellularised human skin tissue transplant: manufacture and properties, E. Rössner, M.D. Smith, B. Petschke, K. Schmidt, M. Vitacolonna, C. Syring, R. von Versen, and P. Hohenberger. Cell and Tissue Banking, Vol, No. 3, pp. 209-217.

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OsteoFleece® is an implantable bio-porcine collagen membrane designed to optimize bone and soft tissue regeneration in oral, implant, and maxillofacial surgery. Its porous structure ensures effective retention of blood, cytokines, and proteases and promotes blood clotting. The parallel-structured fibers distribute the exudate and improve the migration of cells involved in wound healing.

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E

CHARACTERISTICS

- Promotes wound healing
- Supports tissue regeneration
- Completely resorbable
- 100% biocompatible: free of BSE pathogens, bacteria, and viruses
- Anti-inflammatory and analgesic
- Binds a multiple of its own volume in fluid (hemostatic)
- Very well tolerated
- Exhibits very low immunogenicity due to high homology with human collagen

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E

AVAILABLE SIZES

- 20 x 20 mm / 10 pieces per package, Art.Nr.: OF-2020
- The fleece can be cut to your desired sizes.

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RELEVANCE IN EXTRACTION MANAGEMENT

S3-Guideline – Dental surgery in patients on oral anticoagulants.

In the field of dental surgery, postoperative bleeding in particular can pose a serious problem. For patients on oral anticoagulants or antiplatelet agents, there is considerable uncertainty regarding the preparation, performance, and follow-up care associated with dental surgical procedures.

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E

AREAS OF APPLICATION | INDICATIONS

- Promotion of wound healing and support for tissue regeneration
- Hemostatic wound care of the extraction socket
- Reduction of bleeding and post-operative bleeding
- Tissue regeneration of the extraction socket (Socket Preservation)
- Prevention of wound infections following tooth extraction
- Tissue regeneration in internal sinus lifts
- Biopsy sites
- Minor bone defects
- Soft-tissue dehiscences
- Oral wounds

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E OsteoGel consists of non-crosslinked hyaluronic acid and has been specifically designed to treat inflammation of the oral mucosa and to promote or guide wound healing following surgical procedures (GTR). OsteoGel is perfectly suited as a carrier material for granules of the OsteoGraft product portfolio. The low viscosity of the hyaluronic acid liquid gel allows for easy clinical application.

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E **IMPLANTOLOGIE:**

- Activating bone regeneration material with OsteoGel improves its workability and accelerates wound healing.
- Increased osteoblast formation (as well as osteogenesis) to accelerate bone regeneration.
- Supports and accelerates wound healing through its anti-inflammatory properties (inhibition of the production of pro-inflammatory cytokines, e.g., TNF-α).
- Reduced scarring in areas where aesthetics are a priority (the HA-matrix prevents collagen deposition, resulting in less scarring).

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Osteogel hyaluronic acid liquid gel is locally absorbed through the oral mucosa within 6 to 12 hours. To achieve optimal treatment results, it is important that the hyaluronic acid liquid gel is applied to the wound and that the wound is closed with a suture. **E**

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GENERAL ORAL SURGICAL AND DENTAL APPLICATIONS **E**

- Acceleration of healing after vestibuloplasty.
- Reduction of scar tissue.
- Complement to surgical wound closure.
- General inflammatory conditions of the oral mucosa.
- Safety and prevention following surgical treatment through antiseptic and bacteriostatic effects in the wound area.
- Increased fibroblast proliferation.

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BONE GRAFTS BY DIZG

READY TO USE IN 3 STEPS

The conserved grafts are ready to use in a few easy steps.
 All grafts are double sterile packed.
 Allogeneic bone grafts are fully miscible with autologous bone and / or patient blood.

- 1** Take Transplant off from the innermost, **sterile packaging** and place in a **sterile container**.
- 2** **30 minutes** in one **Leave for 30 minutes** for rehydration in an appropriate **physiological medium** (eg in isotonic infusion solution).
 Frozen grafts: Thaw in the same way for 30 minutes.
- 3** **Graft Remove Graft.**
 It is ready now for use.

STRUCTURE GRAFTS

GRANULES AND CHIPS

SOFT TISSUE
 (e.g. Epiflex®)

Please note the expiry date stated on the package and refer to the fact sheet for use and specification on the corresponding transplant.

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■ ARGON Dental Vertriebs GmbH & Co. KG